

Product Datasheet

Bovine Postacrosomal Sheath WW Domain Binding Protein (PAWP) ELISA Kit (orb1664263)

Description	Bovine Postacrosomal Sheath WW Domain Binding Protein (PAWP) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.059 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.