

Product Datasheet

C8B antibody (orb166429)

Description	Rabbit polyclonal antibody to C8B (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	C8B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	67047
Uniprot ID	P07358
NCBI	NP_000057.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

