

Product Datasheet

GRIN2C antibody (orb166437)

Description	Rabbit polyclonal antibody to GRIN2C
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GRIN2C
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	134209
Uniprot ID	Q14957
NCBI	NP_000826.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

