

Product Datasheet

Lrmp antibody (orb166443)

Description	Rabbit polyclonal antibody to Lrmp (C-term).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Irag2
Form/Appearance	Supplied in PBS with 0.09% (W/V) sodium azide.
Concentration	0.44 mg/ml

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal


MW 59588 Da

Uniprot ID [Q60664](#)

NCBI [NP_001268909.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.



Western blot analysis of lysate from mouse neuro-2a cells.