

Product Datasheet

Rat Nidogen 2 (NID2) ELISA Kit (orb1664523)

Description	Rat Nidogen 2 (NID2) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 ng/mL. Sensitivity: 32 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.