

Product Datasheet

Bovine Nerve Growth Factor (NGF) ELISA Kit (orb1664534)

Description	Bovine Nerve Growth Factor (NGF) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 125-8000 pg/mL. Sensitivity: 30 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.