

Product Datasheet

PRSS36 antibody (orb166457)

Description	Rabbit polyclonal antibody to PRSS36
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PRSS36
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	91955
Uniprot ID	Q5K4E3
NCBI	NP_775773.2 , NP_001245219.1 , NP_001245220.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

