

Product Datasheet

Human N-alpha-acetyltransferase 10 (NAA10) ELISA Kit (orb1664613)

Description Human N-alpha-acetyltransferase
10 (NAA10) ELISA Kit

Reactivity Human

Tested ELISA

Applications

Note For research use only

Application notes Detect Range: 31.2-2000 pg/mL.
Sensitivity: 9.64 pg/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
serum, plasma, tissue
homogenates, cell lysates, cell
culture supernates or other
biological fluids.

Expiration Date Please enquire.