

Product Datasheet

Bovine Matrix Metalloproteinase 7 (MMP7) ELISA Kit (orb1664772)

Description	Bovine Matrix Metalloproteinase 7 (MMP7) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.089 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.