

Product Datasheet

HBE1 antibody (orb166478)

Description	Rabbit polyclonal antibody to HBE1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	HBE1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16203
Uniprot ID	P02100
NCBI	NP_005321.1
Dilution Range	WB: 1:8000
Expiration Date	12 months from date of receipt.

