

Product Datasheet

Chicken Matrix Metalloproteinase 13 (MMP13) ELISA Kit (orb1664802)

Description	Chicken Matrix Metalloproteinase 13 (MMP13) ELISA Kit
Reactivity	Gallus
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.63 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.