

Product Datasheet

OR5AC2 antibody (orb166484)

Description	Rabbit polyclonal antibody to OR5AC2 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	OR5AC2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35304
Uniprot ID	Q9NZP5
NCBI	NP_473447.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

