

Product Datasheet

UGT2B28 antibody (orb166496)

Description	Rabbit polyclonal antibody to UGT2B28
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	UGT2B28
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60906
Uniprot ID	Q9BY64
NCBI	NP_001193933.1 , NP_444267.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

