

Product Datasheet

Mouse Leucine Rich Repeat Containing Protein 3C (LRRC3C) ELISA Kit (orb1664988)

Description Mouse Leucine Rich Repeat
Containing Protein 3C (LRRC3C)
ELISA Kit

Reactivity Mouse

**Tested
Applications** ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL.
Sensitivity: 0.068 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
tissue homogenates or other
biological fluids.

Expiration Date Please enquire.