

Product Datasheet

Human Latrophilin 3 (LPHN3) ELISA Kit (orb1665010)

Description Human Latrophilin 3 (LPHN3) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL.
Sensitivity: 0.114 ng/mL. Assay Time:
1-4.5 h. Detection wavelength: 450
nm. Sample Volume: 100 µL. Sample
Type: tissue homogenates, cell lysates
or other biological fluids.

Expiration Date Please enquire.