

Product Datasheet

Monkey Luteinizing Hormone (LH) ELISA Kit (orb1665062)

Description	Monkey Luteinizing Hormone (LH) ELISA Kit
Reactivity	Monkey
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.468-30 ng/mL. Sensitivity: 0.172 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.