

Product Datasheet

Human B-Lymphocyte Activation Antigen B7-2 (LAB7-2) ELISA Kit (orb1665135)

Description Human B-Lymphocyte Activation
Antigen B7-2 (LAB7-2) ELISA Kit

Reactivity Human

Tested ELISA

Applications

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL.
Sensitivity: 0.062 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
tissue homogenates, cell lysates,
cell culture supernates or other
biological fluids.

Expiration Date Please enquire.