

Product Datasheet

Human Kinesin Family, Member 18A (KIF18A) ELISA Kit (orb1665213)

Description Human Kinesin Family, Member 18A (KIF18A) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL. Sensitivity: 0.109 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.

Expiration Date Please enquire.