

Product Datasheet

Human Guanine nucleotide exchange factor VAV3 ELISA Kit (orb1665250)

Description	Human Guanine nucleotide exchange factor VAV3 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: less than 0.078 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 T. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.