

Product Datasheet

4E-BP2 antibody (orb166535)

Description	Rabbit polyclonal antibody against 4E-BP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	EIF4EBP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	12939
Uniprot ID	Q13542
NCBI	NP_004087.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

