

Product Datasheet

Rabbit Interleukin 8 Receptor Alpha (IL8Ra) ELISA Kit (orb1665371)

Description	Rabbit Interleukin 8 Receptor Alpha (IL8Ra) ELISA Kit
Reactivity	Rabbit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.061 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.