

Product Datasheet

ATP6AP2 antibody (orb166540)

Description	Rabbit polyclonal antibody to ATP6AP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ATP6AP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39008
Uniprot ID	O75787
NCBI	NP_005756.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

