

Product Datasheet

Monkey Interleukin 4 (IL4) ELISA Kit (orb1665407)

Description	Monkey Interleukin 4 (IL4) ELISA Kit
Reactivity	Monkey
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 pg/mL. Sensitivity: 1.26 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.