

Product Datasheet

Human Interleukin 4 Induced Protein 1 (IL4I1) ELISA Kit (orb1665412)

Description	Human Interleukin 4 Induced Protein 1 (IL4I1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.052 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.