

Product Datasheet

Rat Alanine aminotransferase 1 ELISA Kit (orb1665433)

Description	Rat Alanine aminotransferase 1 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.78-50 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 T. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.