

Product Datasheet

Guinea pig Interleukin 2 (IL2) ELISA Kit (orb1665440)

Description Guinea pig Interleukin 2 (IL2) ELISA Kit

Reactivity Guinea pig

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 15.6-1000 pg/mL.
Sensitivity: 6.9 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.

Expiration Date Please enquire.