

## Product Datasheet

# Rat Interleukin 23 Receptor (IL23R) ELISA Kit (orb1665460)

|                            |                                                                                                                                                                                                                            |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rat Interleukin 23 Receptor (IL23R)<br>ELISA Kit                                                                                                                                                                           |
| <b>Reactivity</b>          | Rat                                                                                                                                                                                                                        |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                      |
| <b>Note</b>                | For research use only                                                                                                                                                                                                      |
| <b>Application notes</b>   | Detect Range: 0.156-10 ng/mL.<br>Sensitivity: 0.042 ng/mL. Assay Time:<br>1-4.5 h. Detection wavelength: 450<br>nm. Sample Volume: 100 µL. Sample<br>Type: tissue homogenates, cell<br>lysates or other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                                            |