

Product Datasheet

PDCD6IP antibody (orb166549)

Description	Rabbit polyclonal antibody to PDCD6IP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PDCD6IP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	96023
Uniprot ID	Q8WUM4
NCBI	NP_001155901.1 , NP_037506.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

