

Product Datasheet

Monkey Interleukin 1 Receptor Type II (IL1R2) ELISA Kit (orb1665491)

Description	Monkey Interleukin 1 Receptor Type II (IL1R2) ELISA Kit
Reactivity	Monkey
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.213 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.