

Product Datasheet

Rat Interleukin 18 Binding Protein (IL18BP) ELISA Kit (orb1665537)

Description	Rat Interleukin 18 Binding Protein (IL18BP) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.215 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.