

Product Datasheet

Mouse Interleukin 13 Receptor Alpha 2 (IL13Ra2) ELISA Kit (orb1665568)

Description	Mouse Interleukin 13 Receptor Alpha 2 (IL13Ra2) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.126 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.