

Product Datasheet

Human Interleukin 13 Receptor Alpha 2 (IL13Ra2) ELISA Kit (orb1665569)

Description Human Interleukin 13 Receptor
Alpha 2 (IL13Ra2) ELISA Kit

Reactivity Human

Tested ELISA

Applications

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL.
Sensitivity: 0.129 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
serum, plasma, tissue
homogenates, cell lysates, cell
culture supernates or other
biological fluids.

Expiration Date Please enquire.