

## Product Datasheet

# Human Interleukin 12 Receptor Beta 2 (IL12Rb2) ELISA Kit (orb1665578)

|                                |                                                                                                                                                                                                                                   |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>             | Human Interleukin 12 Receptor<br>Beta 2 (IL12Rb2) ELISA Kit                                                                                                                                                                       |
| <b>Reactivity</b>              | Human                                                                                                                                                                                                                             |
| <b>Tested<br/>Applications</b> | ELISA                                                                                                                                                                                                                             |
| <b>Note</b>                    | For research use only                                                                                                                                                                                                             |
| <b>Application notes</b>       | Detect Range: 0.156-10 ng/mL.<br>Sensitivity: 0.055 ng/mL. Assay<br>Time: 1-4.5 h. Detection<br>wavelength: 450 nm. Sample<br>Volume: 100 µL. Sample Type:<br>serum, plasma, tissue<br>homogenates or other biological<br>fluids. |
| <b>Expiration Date</b>         | Please enquire.                                                                                                                                                                                                                   |