



Product Datasheet

Human Activin A / INHBA Protein (orb16656)

Description	Human Activin A, premium grade is expressed from human 293 cells (HEK293). It contains AA Gly 311 - ...
Reactivity	Human
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	Activin A
Preservatives	0.056% TFA in 30% ACN
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 13.0 kDa. As a result of glycosylation, the protein migrates as 14-15 kDa under reducing (R) condition, and 23-25 kDa under non-reducing (NR) condition (SDS-PAGE).
Purity	95%
Source	Human Activin A, premium grade (orb16656) is expressed from human 293 cells (HEK293). It contains AA Gly 311 - Ser 426 (Accession # AAH07858.1).
MW	13.0 kDa
Uniprot ID	P08476
NCBI	AAH07858.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Activin A pro...



Immobilized Human Activin A, premium gra...