

Product Datasheet

Rat Interleukin 10 Receptor Beta (IL10Rb) ELISA Kit (orb1665600)

Description	Rat Interleukin 10 Receptor Beta (IL10Rb) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.065 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.