

Product Datasheet

Mouse Insulin Like Growth Factor 1 (IGF1) ELISA Kit (orb1665676)

Description	Mouse Insulin Like Growth Factor 1 (IGF1) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 29 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.