



Product Datasheet

CLDN14 antibody (orb166573)

Description	Rabbit polyclonal antibody to CLDN14 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CLDN14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25699
Uniprot ID	O95500
NCBI	NP_652763.1, NP_001139549.1, NP_001139550.1, NP_036262.1, NP_001139551.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis of human kidney and liver ti...