

Product Datasheet

MBNL1 antibody (orb166577)

Description	Rabbit polyclonal antibody to MBNL1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MBNL1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41817
Uniprot ID	Q9NR56
NCBI	NP_997176.1 , NP_997180.1 , NP_997175.1 , NP_066368.2 , NP_997177.1 , NP_997179.1 , NP_997178.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

