

Product Datasheet

HNRNPH1 antibody (orb166589)

Description	Rabbit polyclonal antibody to HNRNPH1 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	HNRNPH1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49229
Uniprot ID	P31943
NCBI	NP_005511.1, NP_001244222.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

