

Product Datasheet

Human Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 2C (GRIN2C) ELISA Kit (orb1666010)

Description Human Glutamate Receptor,
Ionotropic, N-Methyl-D-Aspartate
2C (GRIN2C) ELISA Kit

Reactivity Human

**Tested
Applications** ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL.
Sensitivity: 0.118 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
tissue homogenates or other
biological fluids.

Expiration Date Please enquire.