

## Product Datasheet

# Mouse Glycoprotein Ib Beta Polypeptide, Platelet (GP1Bb) ELISA Kit (orb1666066)

|                            |                                                                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Glycoprotein Ib Beta Polypeptide, Platelet (GP1Bb) ELISA Kit                                                                                     |
| <b>Reactivity</b>          | Mouse                                                                                                                                                  |
| <b>Tested Applications</b> | ELISA                                                                                                                                                  |
| <b>Note</b>                | For research use only                                                                                                                                  |
| <b>Application notes</b>   | Detect Range: 31.2-2000 pg/mL. Sensitivity: 14.7 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: plasma. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                        |