

Product Datasheet

Txnip antibody (orb166623)

Description	Rabbit polyclonal antibody to Txnip
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	Txnip
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	44363
Uniprot ID	Q8BG60
NCBI	NP_076208.2 , NP_001009935.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

