

Product Datasheet

Rat Fibrinogen Gamma (FGg) ELISA Kit (orb1666364)

Description	Rat Fibrinogen Gamma (FGg) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 ng/mL. Sensitivity: 27 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: plasma.
Expiration Date	Please enquire.