

Product Datasheet

Bovine Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit (orb1666393)

Description Bovine Fibroblast Growth Factor 2,
Basic (FGF2) ELISA Kit

Reactivity Bovine

Tested ELISA

Applications

Note For research use only

Application notes Detect Range: 15.6-1000 pg/mL.
Sensitivity: 6.5 pg/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
serum, plasma, tissue
homogenates, cell lysates, cell
culture supernates or other
biological fluids.

Expiration Date Please enquire.