

Product Datasheet

Rat Fibroblast Growth Factor 11 (FGF11) ELISA Kit (orb1666417)

Description	Rat Fibroblast Growth Factor 11 (FGF11) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.2 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.