

Product Datasheet

Porcine Factor Related Apoptosis Ligand (FASL) ELISA Kit (orb1666478)

Description	Porcine Factor Related Apoptosis Ligand (FASL) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 23.4-1500 pg/mL. Sensitivity: 5.2 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.