

Product Datasheet

Human Fatty Acid Binding Protein 2, Intestinal (FABP2) ELISA Kit (orb1666516)

Description	Human Fatty Acid Binding Protein 2, Intestinal (FABP2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.057 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.