

Product Datasheet

Mouse Coagulation Factor IX (F9) ELISA Kit (orb1666520)

Description	Mouse Coagulation Factor IX (F9) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.067 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: plasma.
Expiration Date	Please enquire.