

Product Datasheet

PPP2CA antibody (orb166654)

Description	Rabbit polyclonal antibody to PPP2CA
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PPP2CA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35594
Uniprot ID	P67775
NCBI	NP_002706.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

