

Product Datasheet

Human Enhancer Of Zeste Homolog 1 (EZH1) ELISA Kit (orb1666540)

Description	Human Enhancer Of Zeste Homolog 1 (EZH1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.113 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.